

Product datasheet for **RC218237A1V**

Human PKI alpha (PKIA) (NM_006823) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PKI alpha (PKIA) (NM_006823) AAV Particle
Tag: Myc-DDK
Symbol: PKI alpha
Synonyms: PRKACN1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC218237 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACTGATGTGGAACTACATATGCAGATTTATTGCTTCAGGAAGAACAGGTAGAAGAAATGCAATAC
ATGATATCCTGGTTTCTCTGCAAGTGGCAACAGCAATGAATTAGCCTTGAAATTAGCAGGTCTTGATAT
CAACAAGACAGAAGGTGAAGAAGATGCACAACGAAGTTCTACAGAACAAAGTGGGAAGCCAGGGAGAA
GCAGCAAAATCTGAAAGC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC218237 protein sequence
Red=Cloning site Green=Tags(s)

MTDVETTYADFIASGRTGRRNAIHDILVSSASGNSNELALKLAGLDINKTEGEEDAQRSSTEQSGEAQGE
AAKSES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_006823
ORF Size: 228 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006823.2</u>
RefSeq Size:	4215 bp
RefSeq ORF:	231 bp
Locus ID:	5569
UniProt ID:	<u>P61925</u>
Cytogenetics:	8q21.13
MW:	8 kDa