

Product datasheet for **RC218251A1V**

Human Sterol carrier protein 2 (SCP2) (NM_001007100) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Sterol carrier protein 2 (SCP2) (NM_001007100) AAV Particle
Tag:	Myc-DDK
Symbol:	Sterol carrier protein 2
Synonyms:	NLTP; NSL-TP; SCOX; SCP-2; SCP-CHI; SCP-X; SCPX
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC218251 representing NM_001007100 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC ATGGGTTTTCCGAAGCCGCCAGAACTCATCAAATTGAAGCTGTTCCAACCAGCTCTGCAAGTGATGGATTTAAGGCAAATCTTGTTTTAAGGAGATTGAGAAGAACTTGAAGAGGAAGGGGAACAGTTTGTGAAGAAATCGGTGGTATTTTGCCTTCAAGGTGAAAGATGGCCCTGGGGGTAAAGAGGCCACCTGGTGGTGGATGTGAAGAATGGCAAAGGATCAGTGCTTCTTAAGTCAAGTGAAGAAGGCTGACTGCACAATCACAATGGCTGACTCAGACTTCTGGCTTTAATGACTGGTAAATGAATCCTCAGTCGGCCTTCTTCAAGGCAAATTGAAATCACTGGCAACATGGGTCTCGCTATGAAGTTACAAAATCTTCAGCTTCAGCCAGGCAACGCTAAGCTC ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC218251 representing NM_001007100 Red=Cloning site Green=Tags(s) MGFPEAARTHQIEAVPTSSASDGFKANLVFKEIEKKLEEEGEQFVKKIGGIFAFKVKDGPGGKEATWVVDVKNKGKSVLPNSDKKADCTITMADSDFLALMTGKMNPQSAFFQGKLKITGNMGLAMKLQNLQLQPGNAKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001007100


[View online »](#)

ORF Size:	420 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001007100.1</u>
RefSeq Size:	1438 bp
RefSeq ORF:	423 bp
Locus ID:	6342
UniProt ID:	<u>P22307</u>
Cytogenetics:	1p32.3
MW:	15.08 kDa