

## Product datasheet for **RC218270L2V**

### SCARB1 (NM\_001082959) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SCARB1 (NM_001082959) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCARB1
Synonyms:	CD36L1; CLA-1; CLA1; HDLQTL6; SR-BI; SRB1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001082959
ORF Size:	1518 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218270).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001082959.1</a> , <a href="#">NP_001076428.1</a>
RefSeq Size:	2630 bp
RefSeq ORF:	1521 bp
Locus ID:	949
UniProt ID:	<a href="#">Q8WTV0</a>
Cytogenetics:	12q24.31
Protein Families:	Druggable Genome, Transmembrane
MW:	56 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a plasma membrane receptor for high density lipoprotein cholesterol (HDL). The encoded protein mediates cholesterol transfer to and from HDL. In addition, this protein is a receptor for hepatitis C virus glycoprotein E2. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2019]