

Product datasheet for **RC218288**

C5ORF44 (TRAPPC13) (NM_001093756) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C5ORF44 (TRAPPC13) (NM_001093756) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C5ORF44
Synonyms:	C5orf44; SHLD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC218288 representing NM_001093756 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGAAGTGAATCCCCCTAAACAGGAGCACCTGCTGGCGCTAAAAGTGATGCGGCTGACTAAGCCTACTT
TATTCACCAATATCCAGTAACATGTGAAGAGAAAGACTTACCTGGAGATCTCTTAAACCAGCTGATGAG
AGATGATCCTTCAACCGTTAATGGTGCAGAAGTTTAAATGTTGGGAGAAATGCTGACTTTACCACAGAAT
TTTGGGAATATATTTTGGGAGAGACCTTTCCAGTTATATCAGCGTTCATAATGATAGCAATCAAGTTG
TAAAAGACATATTAGTAAAAGCTGATCTTCAGACAAGTTCTCAGCGTTTAAATCTTTCAGCCTCCAATGC
TGCACTGGCTGAACTTAAACCGGATTGTTGTATTGATGATGTCATACATCATGAAGTCAAAGAAATTGGA
ACACACATCTTGGTATGTGCTGTGAGTTATACTCACTCAGGCTGGAGAAAAATGTATTTAGAAAAATCT
TCAAAATTCAGGTTCTCAAACCATTTGGATGTGAAAACCAATTTTACAATGCAGAGACTGATGAAGTATT
TCTGGAAGCCAGATTGAGAATATGACAACCTCACCTATGTTTATGGAGAAGGTTTCACTGGAGCCATCT
ATTATGTACAATGTAAACAGAATTAATTCAGTCAGCCAAGCTGGAGAATGTGTGTCTACGTTTGGGTCAA
GAGCATATTTGCAACCAATGGATACACGCCAGTACTTATACTGCCTAAAGCCAAAGAATGAATTTGCAGA
AAAAGCAGGCATCATTAAAGGGAGTAACAGTAATTGGAAAATTGGATATAGTATGGAAAACAAATCTAGGT
GAAAGGGGAAGGTTACAGACCAGCCAACCTCAAAGAATGGCTCCAGGTTATGGAGATGTTAGGTTGTCTT
TGGAGGCAATACCAGATACCGTAAACCTTGAAGAACCTTTTCATATTACCTGTAAAATAACAACTGCAG
TGAAAGGACTATGGATCTGGTTTTGAAAATGTGAATACCAATCCATCCACTGGTGTGGAATTTACAGGA
AGACAGCTGGGAAAGCTGCATCCAAGTTCTCGCTCTGTCTTGCCTTACTCTGCTTTCTTCAGTACAGG
GACTGCAAAGCATCTCTGGCTTAAGACTAACAGACACATTCTTAAAGAGAACATATGAATATGATGACAT
CGCACAGTCTGTGTGGTATCTTCTGCCATTAAAGTGAAAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC218288 representing NM_001093756
 Red=Cloning site Green=Tags(s)


MEVNPPKQEHLLALKVMRLTKPTLFTNIPVTCEEKDLPGDLFNQLMRDDPSTVNGAEVLMGEMLTLPQN
 FGNIFLGETFSSYISVHNSNQVVKDILVKADLQTSSQRLNLSASNAVAELKPDCIDDVIHHEVKEIG
 THILVCAVSYTTQAGEKMYFRKFFKFQVLKPLDVKTKFYNAETDEVFLEAQIQNMTTSPMFMEKVSLEPS
 IMYNVTELNVSQAGECVSTFGSRAYLQPMQTRQYLCLKPKNEFAEKAGIIGVTVIGKLDIVWKTNLG
 ERGRLQTSQLQRMAGPYGDVRLSLEAIPDTVNLEEPFHITCKITNCERTMDLVLEMNTNSIHWCGISG
 RQLGKLHPSSSLCLALTLLSSVQGLQSIGSLRLTDTFLKRTYEYDDIAQVCVSSAIKVES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1304_d10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001093756

ORF Size: 1233 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001093756.1](#), [NP_001087225.1](#)

RefSeq Size: 3099 bp

RefSeq ORF: 1236 bp


Locus ID: 80006

UniProt ID: [A5PLN9](#)


Cytogenetics: 5q12.3

MW: 45.8 kDa


Product images:



Circular map for RC218288



Western blot validation of overexpression lysate (Cat# [LY420689]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC218288 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified TRAPPC13 protein (Cat# [TP318288]). The protein was produced from HEK293T cells transfected with TRAPPC13 cDNA clone (Cat# RC218288) using MegaTran 2.0 (Cat# [TT210002]).