

## Product datasheet for **RC218298L2V**

### **C14orf130 (UBR7) (NM\_175748) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C14orf130 (UBR7) (NM_175748) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C14orf130
Synonyms:	C14orf130; LICAS
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_175748
ORF Size:	1275 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218298).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175748.1</a>
RefSeq Size:	1576 bp
RefSeq ORF:	1278 bp
Locus ID:	55148
UniProt ID:	<a href="#">Q8N806</a>
Cytogenetics:	14q32.12
Protein Families:	Druggable Genome
MW:	47.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a UBR box-containing protein that belongs to the E3 ubiquitin ligase family. The protein also contains a plant homeodomain (PHD) in the C-terminus. In mammals, the encoded protein recognizes N-degrons, the destabilizing N-terminal residues of short-lived proteins, which results in ubiquitinylation, and proteolysis via the proteasome. [provided by RefSeq, Jul 2016]