

## Product datasheet for **RC218309L1V**

### USP1 (NM\_001017416) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	USP1 (NM_001017416) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP1
Synonyms:	UBP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001017416
ORF Size:	2355 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218309).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001017416.1</a> , <a href="#">NP_001017416.1</a>
RefSeq Size:	3424 bp
RefSeq ORF:	2358 bp
Locus ID:	7398
UniProt ID:	<a href="#">O94782</a>
Cytogenetics:	1p31.3
Protein Families:	Protease
MW:	88 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases that is a deubiquitinating enzyme (DUB) with His and Cys domains. This protein is located in the cytoplasm and cleaves the ubiquitin moiety from ubiquitin-fused precursors and ubiquitinated proteins. The protein specifically deubiquitinates a protein in the Fanconi anemia (FA) DNA repair pathway. Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008]