

Product datasheet for **RC218322A1V**

Human DUSP13 (NM_001007272) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DUSP13 (NM_001007272) AAV Particle
Tag:	Myc-DDK
Symbol:	DUSP13
Synonyms:	BEDP; DUSP13A; DUSP13B; MDSP; SKRP4; TMDP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC218322 representing NM_001007272 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTGAGACCTCTCTCCAGAGCTGGGGGAGAGGACAAAGCCACGCCTTGCCCCAGCATCTGGAGC
TGGAGGAGCTCCTGCGGCAGGGAAGTCTTCTGCAGCCGTGTGGACGAAGTTTGGCCCAACCTTTTCAT
AGGAGATGCGATGGACTCACTGCAGAAGCAGGACCTCCGGAGGCCCAAGATCCATGGGGCAGTCCAGGCA
TCTCCCTACCAGCCGCCACATTGGCTTCGCTGCAGCGCTTGTGTGGGTCCGTACGGCTGCCACACTGA
ACCATATCGATGAGGTCTGGCCAGCCTCTTCTGGGAGATGCGTACGCAGCCCGGACAAAGCAAGT
GATCCAGCTGGGAATCACCCACGTTGTGAATGCCGCTGCAGGCAAGTTCCAGGTGGACACAGGTGCCAAA
TTCTACCGTGGAATGTCCTGGAGTACTATGGCATCGAGGCGGACGACAACCCCTTCTTCGACCTCAGTG
TCTACTTTCTGCCTGTTGCTCGATACATCCGAGCTGCCCTCAGTGTTCCCAAGGCCGCGTGCTGGTACA
CTGTGCCATGGGGTAAGCCGCTCTGCCACACTTGCTCTGGCCTTCTCATGATCTGTGAGAACATGACG
CTGGTAGAGGCCATCCAGACGGTGCAGGCCACCGCAATATCTGCCCTAACTCAGGCTTCTCCGGCAGC
TCCAGGTTCTGGACAACCGACTGGGGCGGGAGACGGGGCGTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC218322 representing NM_001007272 Red=Cloning site Green=Tags(s) MAETSLPELGGEDKATPCPSILELEELLRAGKSSCSRVDVWPNFLFIGDAMDSLQKQDLRRPKIHGAVQA SPYQPPTLASLQRLWVRQAATLNHIDEVWPSLFLGDYAARDKSKLIQLGITHVVNAAAGKFQVDTGAK FYRGSLEYGYIEADDNPFDFLSVYFLPVARYIRAALSVPQGRVLVHCAMGVSRSATLVLAFLMICENMT LVEAIQTVQAHNRNICPNSGFLRQLQVLNRLGRETGRF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001007272
ORF Size:	744 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001007272.1</u>
RefSeq Size:	1063 bp
RefSeq ORF:	747 bp
Locus ID:	51207
UniProt ID:	<u>Q6B8I1</u>
Cytogenetics:	10q22.2
MW:	27.3 kDa