

Product datasheet for **RC218337**

C16orf13 (METTL26) (NM_032366) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C16orf13 (METTL26) (NM_032366) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C16orf13
Synonyms:	C16orf13; JFP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC218337 representing NM_032366 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGTGGCGGCGGCCGCGGAGCGGAACAAGGATCCCATCTTGACGTGCTGCGGCAGTACCTGGATC
CGGCCCAGCGTGGCGTCCGCGTCTCGAGGTGGCTCGGGCTCCGGCCAGCACGCAGCGCACTTCGCGCG
GGCCTTCCCCCTGGCCGAGTGGCAGCCGTGGACGTGGACAGCGCTGCCTGGACAGCATCGCGGCCACC
ACGCAAGCCCAGGGCCTGACCAACGTGAAGGCCCGCTACACCTGGACGTGACGTGGGGCTGGGAGCACT
GGGGCGGGATCCTGCCACAGTCGCTGGACCTGTTGCTCTGCATCAACATGGCCCATGTACGCCCCCTGCG
CTGCACGGAGGGGCTCTTCAGAGCAGCAGGACACCTGCTCAAACCCAGGGCCCTGCTCATCACCTACGGG
CCCTATGCCATCAATGGGAAGATCTCCCCCAGAGCAACGTGGACTTTGACCTGATGCTCAGATGCAGGA
ACCCAGAATGGGGGCTTCGGGACACAGCCCTCCTGGAGGACCTGGGAAAGGCCAGTGGCCTGCTCCTGGA
GAGGATGGTGGACATGCCAGCCAACAACAAATGCCTGATCTTCGGAAAAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC218337 representing NM_032366 Red=Cloning site Green=Tags(s)
-------------------	--

MLVAAAERNKDPILHVLRLQYLDPAQRGVRVLEVASGSGQHAHFARAFPLAEWQPSDQVDRCLDSIAAT
TQAQGLTNVKAPLHLDVTWGEHWGGILPQSLDLLLCINMAHVSPLRCTEGLFRAAGHLLKPRALLITYG
PYAINGKISPQSNVDFDMLRCRNPWGLRDTALLEDLGKASGLLLERMVDMPPANNKCLIFRKN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mk6707_a06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_032366

ORF Size: 612 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_032366.5](#)


RefSeq Size: 854 bp

RefSeq ORF: 615 bp
 Locus ID: 84326
 UniProt ID: [Q96S19](#)
 Cytogenetics: 16p13.3
 MW: 22.4 kDa


Product images:



Circular map for RC218337



Western blot validation of overexpression lysate (Cat# [LY410160]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC218337 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified METTL26 protein (Cat# [TP318337]). The protein was produced from HEK293T cells transfected with METTL26 cDNA clone (Cat# RC218337) using MegaTran 2.0 (Cat# [TT210002]).