

Product datasheet for **RC218337A1V**

Human **C16orf13 (METTL26) (NM_032366) AAV Particle**

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: C16orf13

Synonyms: C16orf13; JFP2

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC218337 representing NM_032366
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGCTGGTGGCGCGCGCGGAGCGGAACAAGGATCCCATCTTGACGTGCTGCGGCAGTACCTGGATC
CGGCCAGCGTGGCGTCCGCTCCTCGAGGTGGCTCGGGCTCCGGCCAGCACGCAGCGCACTTCGCGCG
GGCCTTCCCCCTGGCCGAGTGGCAGCCGTGGGACGTGGACAGCGCTGCCTGGACAGCATCGCGGCCACC
ACGCAAGCCAGGGCCTGACCAACGTGAAGGCCCGCTACACCTGGACGTGACGTGGGGCTGGGAGCACT
GGGGCGGGATCCTGCCACAGTCGCTGGACCTGTTGCTCTGCATCAACATGGCCCATGTCAGCCCCCTGCG
CTGCACGGAGGGGCTCTCAGAGCAGCAGGACCTGCTCAAACCCAGGGCCCTGCTCATCACCTACGGG
CCCTATGCCATCAATGGGAAGATCTCCCCCAGAGCAACGTGGACTTTGACCTGATGCTCAGATGCAGGA
ACCCAGAATGGGGGCTTCGGGACACAGCCCTCCTGGAGGACCTGGGAAAGGCCAGTGGCCTGCTCCTGGA
GAGGATGGTGGACATGCCAGCCAACAACAAATGCCTGATCTTCGGAAAAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC218337 representing NM_032366
Red=Cloning site Green=Tags(s)

MLVAAAERNKDPILHVLRYQLDPAQRGVRVLEVASGSGQHAHFARAFPLAEWQPSDQDQRC LDSIAAT
TQAQGLTNVKAPLHLDVTWGEHWGGILPQSLDLLLCINMAHVSPLRCTEGLFRAAGHLLKPRALLITYG
PYAINGKISPQSNVDFDLMLRCRNP EWGLRDTALLEDLGKASGLLLERMVDMPPANNKCLIFRKN

TRTRPLE**Q**KLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_032366
ORF Size:	612 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_032366.3
RefSeq Size:	854 bp
RefSeq ORF:	615 bp
Locus ID:	84326
UniProt ID:	Q96S19
Cytogenetics:	16p13.3
MW:	22.4 kDa