

## Product datasheet for **RC218339A1V**

### Human OR4F6 (NM\_001005326) AAV Particle

#### Product data:

|                           |                                                                                |
|---------------------------|--------------------------------------------------------------------------------|
| Product Type:             | AAV Particles                                                                  |
| Product Name:             | Human OR4F6 (NM_001005326) AAV Particle                                        |
| Tag:                      | Myc-DDK                                                                        |
| Symbol:                   | OR4F6                                                                          |
| Synonyms:                 | OR4F12; OR15-15                                                                |
| Mammalian Cell Selection: | None                                                                           |
| Vector:                   | pAAV-AC-Myc-DDK (PS100089)                                                     |
| ORF Nucleotide Sequence:  | >RC218339 representing NM_001005326<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGATGAAGCCAATCACTCTGTGGTCTCTGAGTTTGTGTTCTGGGACTCTCTGACTCGCGGAAGATCC  
AGCTCCTCCTCTCTCTCTTTTCTCAGTGTTCTATGTGTCAAGCCTGATGGGAAATCTCCTCATTGTGCT  
AACTGTGACCTCTGACCCTCGTTTACAGTCCCCATGTACTTCCTGCTGGCCAACCTTTCCATCATCAAT  
TTGGTATTTTGTCTCCACAGCTCCCAAGATGATTTATGACCTTTTCAGGAAGCACAAGACCATCTCTT  
TTGGGGGCTGTGTAGTTCAGATCTTCTTTATCCATGCAGTTGGGGGAAGTGAAGTGGTGTCTCATAGC  
CATGGCTTTTGACCGATATGTGGCCATATGTAAGCCTCTCCACTACCTGACCATCATGAACCCACAAGG  
TGCATTTTGTCTTTAGTCATTTCTGGATTATAGGTATTATCACTCAGTGATTGAGTTGGCTTTTGTG  
TAGACCTGCTGTTCTGTGGCCCTAATGAATTAGATAGTTTCTTTTGTGATCTTCCTCGATTTATCAAAC  
GGCTTGCATAGAGACCTACACATTGGGATTCATGGTACTGCCAATAGTGGATTATTTCTCTGGCTTCT  
TTTTTAATTCATAATCTCTTACATCTTTATTTTGGTGACTGTTGAGAAAAATCTTCAGGTGGTATAT  
TCAAGGCTTTCTCTATGCTGTCAGCTCATGTCTTGTGGTGGTTTGGTCTTTGGGCCATTAATCTTTT  
CTATATTTTCCATTTCCACATCACATCTTGATAAATTCCTTGCCATCTTTGATGCAGTTATCACTCCC  
GTTTTGAATCCAGTCATCTATACTTTTAGAAATAAGAGATGATGGTGGCAATGAGAAGACGATGCTCTC  
AGTTTGTGAATTACAGTAAATCTTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Protein Sequence:</b> | <p>&gt;RC218339 representing NM_001005326</p> <p>Red=Cloning site Green=Tags(s)</p> <p>MDEANHSVVSEFVFLGLSDSRKIQLLLFLFFSVFYVSSLMGNLLIVLTVTSDPRLQSPMYFLLANLSIIN<br/> LVFCSSTAPKMIYDLFRKHKTISFGGCVVQIFFIHAVGGTEMVLLIAMAFDRYVAICKPLHYLTIMNPQR<br/> CILFLVISWIIIGIHSVIQLAFVVDLLFCGPNELDSFFCDLPRFIKLACIETYTLGFMVTANSGFISLAS<br/> FLILIIISYIFILVTVQKKSSGGIFKAFSMLSAHVIVVVLVFGPLIFFYIFPFPTSHLDKFLAIFDAVITP<br/> VLNPVIYTFRNKEMMVAMRRRCSQFVNYSKIF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p> |
| <b>Species:</b>          | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Serotype:</b>         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>ACCN:</b>             | NM_001005326                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>ORF Size:</b>         | 936 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Buffer:</b>           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Stability:</b>        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.         |
| <b>RefSeq:</b>           | <u>NM_001005326.1, NP_001005326.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>RefSeq Size:</b>      | 939 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>RefSeq ORF:</b>       | 939 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Locus ID:</b>         | 390648                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>UniProt ID:</b>       | <u>Q8NGB9</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Cytogenetics:</b>     | 15q26.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>MW:</b>               | 35.2 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |