

Product datasheet for **RC218351A1V**

Human ZNF137 (ZNF137P) (NM_003438) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: ZNF137

Synonyms: pHZ-30, MGC119990, MGC119991

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC218351 representing NM_003438
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAATGTGGCAAGGTTTTAGTAGAAAAGCACACCTTGCATGTCATCATAGACTTCATACTGTCTAAGG
TTTCTAATCAACAATCAAACCTTGCACAACATCAGAGAGTTTATACTGGAGAGAAACCTTACAAGTGTA
TGAGTGGGGCAAAGCCTTAAGTGGGAAGTCGTCACCTTTTTATCATCAAGCAATCCATGGTGTAGGGAAA
CTTTGCAATGTAATGATTGTCAAAAGTCTTCAGTAATGCTACAACCATTGCAATCACTGGAGAAATCC
ATAATGAAGACAGATCTTACAAGTGAATAAATGTGGTAAATTTTCAGACATCGATCATATCTTGCACT
TTATCAGCGAACTCATACTGGAGAGAAACCTTACAAATATCATGACTGTGGCAAGGTCTTCAGTCAAGCT
TCATCCTATGCAAAACATAGGAGAATTCATACAGGAGAGAAACCTCACAAAGTGTGATGATTGTGGCAAAG
TCTTGACTTCACGTTACACCTCATTAGACATCAGAGAATCCATACTGGACAGAAATCTTACAAATGTCT
TAAGTGTGGCAAGGTCTTCAGTCTGTGGCACTCCATGCAGAACATCAGAAAATTCATTTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC218351 representing NM_003438
Red=Cloning site Green=Tags(s)

MNVARFLVEKHTLHVIIIDFILSKVSNQQSNLAQHQRVYTGEKPYKCNWKGKALSGKSSLFYHQAIHGVGK
LCKCNDCHKVFSNATTIANHWRIHNEEDRSYKCNKCGKIFRHRSYLAVYQRTHTGEKPYKYHDCGKVFSSQA
SSYAKHRRIHTEKPKCDDCGKVLTSRSHLIRHQRIHTGQKSYKCLKCGKVFSLWALHAHQKIHF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_003438
ORF Size:	621 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003438.1
RefSeq Size:	2553 bp
RefSeq ORF:	623 bp
Locus ID:	7696
Cytogenetics:	19q13.41
MW:	23.9 kDa