

Product datasheet for **RC218356L4V**

ZNF750 (NM_024702) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF750 (NM_024702) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF750
Synonyms:	ZFP750
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024702
ORF Size:	2169 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218356).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024702.1
RefSeq Size:	3218 bp
RefSeq ORF:	2172 bp
Locus ID:	79755
UniProt ID:	Q32MQ0
Cytogenetics:	17q25.3
Domains:	zf-C2H2
MW:	77.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein with a nuclear localization site and a C2H2 zinc finger domain. Mutations in this gene have been associated with seborrhea-like dermatitis with psoriasiform elements. [provided by RefSeq, Jul 2008]