

Product datasheet for RC218357A1V

Human TEX261 (NM_144582) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: TEX261

Synonyms: TEG-261

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC218357 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGTGGTTCATGTACCTGCTGAGCTGGCTGTCGCTCTTCATCCAGGTGGCCTTCATCAGCTGGCTGTCG
CGGCTGGACTCTATTACCTGGCAGAACTGATAGAAGAATACACAGTGGCCACCAGCAGGATCATAAATA
CATGATCTGGTCTCCACCGCTGTACTGATTGGCCTCTACGTCTTTGAGCGCTTCCCCACCAGCATGATT
GGAGTGGGCTATTCACCAACCTCGTCTACTTTGGCCTCTCCAGACCTTCCCCTTCATCATGCTGACCT
CGCCTAACTTCATCCTGTCGTGTGGACTAGTGGTGGTGAATCATTACCTAGCATTTCAGTTTTTGCAGA
AGAATATTATCCCTTCTCAGAGGTCTGGCCTATTTCACTTTCTGCCTGTGGATAATTCCGTTTGCCTTT
TTTGTGTCACCTTTTCGGCCGGGGAGAACGTCCTGCCCTCTACCATGCAGCCAGGAGATGATGCTGCTCCA
ATTATTTACCAAAGGCAAGCGGGGCAACGCTTAGGGATCCTGGTTGCTTCTCCTTCATCAAAGAGGC
CATTCTACCGAGTCGTGAGAAGATATAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC218357 protein sequence
Red=Cloning site Green=Tags(s)

MWFMYLLSWLSLFIQVAFITLAVAAGLYLAELIEEYTVATSRIIKYMIWFSTAVLIGLYVFERFPTSMI
GVGLFTNLVYFGLLQTFPFIIMLTSPNFILSCGLVVVNHLYAFQFFAEYYPFSEVLAYFTFCLWIIPFAF
FVLSAGENVLPSTMQPGDDVVSNYFTKGRGKRLGILVVSFIIKEAILPSRQKIY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_144582
ORF Size:	588 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_144582.2
RefSeq Size:	3381 bp
RefSeq ORF:	591 bp
Locus ID:	113419
UniProt ID:	Q6UWH6
Cytogenetics:	2p13.3
MW:	22.5 kDa