

Product datasheet for RC218368L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GUCY2D (NM_000180) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: GUCY2D (NM_000180) Human Tagged Lenti ORF Clone

Tag: mGFP
Symbol: GUCY2D

Synonyms: CACD1; CG-E; CORD5; CORD6; CSNB1I; CYGD; GUC1A4; GUC2D; LCA; LCA1; RCD2; retGC;

RETGC-1; ROS-GC1; ROSGC

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide


The ORF insert of this clone is exactly the same as(RC218368).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_000180 **ORF Size:** 3309 bp





GUCY2D (NM_000180) Human Tagged Lenti ORF Clone - RC218368L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000180.1</u>

RefSeq Size: 3621 bp
RefSeq ORF: 3312 bp
Locus ID: 3000
UniProt ID: Q02846

Cytogenetics: 17p13.1

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Olfactory transduction, Purine metabolism

MW: 119.9 kDa


Gene Summary: This gene encodes a retina-specific guanylate cyclase, which is a member of the membrane

guanylyl cyclase family. Like other membrane guanylyl cyclases, this enzyme has a hydrophobic amino-terminal signal sequence followed by a large extracellular domain, a single membrane spanning domain, a kinase homology domain, and a guanylyl cyclase catalytic domain. In contrast to other membrane guanylyl cyclases, this enzyme is not activated by natriuretic peptides. Mutations in this gene result in Leber congenital amaurosis

and cone-rod dystrophy-6 diseases. [provided by RefSeq, Dec 2008]



Product images:



Circular map for RC218368L4