

Product datasheet for **RC218377L3V**

Kallikrein 10 (KLK10) (NM_145888) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | Lentiviral Particles |
| Product Name: | Kallikrein 10 (KLK10) (NM_145888) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Kallikrein 10 |
| Synonyms: | NES1; PRSSL1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_145888 |
| ORF Size: | 828 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC218377). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_145888.1 |
| RefSeq Size: | 2985 bp |
| RefSeq ORF: | 831 bp |
| Locus ID: | 5655 |
| UniProt ID: | O43240 |
| Cytogenetics: | 19q13.41 |
| Protein Families: | Druggable Genome, Protease, Secreted Protein |
| MW: | 30.1 kDa |



[View online »](#)

Gene Summary:

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Its encoded protein is secreted and may play a role in suppression of tumorigenesis in breast and prostate cancers. Alternate splicing of this gene results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]