

Product datasheet for **RC218383A1V**

Human MMP28 (NM_001032278) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human MMP28 (NM_001032278) AAV Particle

Tag: Myc-DDK

Symbol: MMP28

Synonyms: EPILYSIN; MM28; MMP-25; MMP-28; MMP25

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC218383 representing NM_001032278
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCGCGCGCGTGGCCTCCTGCTGCGCGCCCTGCAGCTGCTACTGTGGGGCCACCTGGACGCCCAGC
CCGCGGAGCGCGGAGGCCAGGAGCTGCGCAAGGAGCGGAGGCATTCTAGAGAAGTACGGATACCTCAA
TGAACAGGTCCCCAAGCTCCACCTCCACTCGATTGAGCGATGCCATCAGAGCGTTTCAGTGGGTGTCC
CAGCTACCTGTGAGCGGCGTGTGGACCGCGCCACCCTGCGCCAGATGACTCGTCCCCGCTGCGGGGTTA
CAGATACCAACAGTTATGCGGCCTGGGCTGAGAGGATCAGTGACTTGTGCTAGACACCGGACCAAAAT
GAGGCGTAAGAAACGCTTTGCAAGCAAGGTGAGCACTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC218383 representing NM_001032278
Red=Cloning site Green=Tags(s)

MVARVGLLLRALQLLLWGHLDAQPAERGGQELRKEAEAFLEKYGYLNEQVPKAPTSTRFSDAIRAFQWVS
QLPVSGVLDRLTRQMTTRPCGVTDNTSYAAWAERISDLFARHRTKMRRKKRFAKQGEHC

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001032278



[View online »](#)

ORF Size:	390 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001032278.1</u>
RefSeq Size:	1042 bp
RefSeq ORF:	393 bp
Locus ID:	79148
UniProt ID:	<u>Q9H239</u>
Cytogenetics:	17q12
MW:	14.93 kDa