

Product datasheet for **RC218401L4V**

GCIP interacting protein p29 (SYF2) (NM_207170) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GCIP interacting protein p29 (SYF2) (NM_207170) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GCIP interacting protein p29
Synonyms:	CBPIN; fSAP29; NTC31; P29
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207170
ORF Size:	603 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218401).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207170.2
RefSeq Size:	1644 bp
RefSeq ORF:	606 bp
Locus ID:	25949
UniProt ID:	O95926
Cytogenetics:	1p36.11
Protein Pathways:	Spliceosome



[View online »](#)

MW: 23.4 kDa

Gene Summary: This gene encodes a nuclear protein that interacts with cyclin D-type binding-protein 1, which is thought to be a cell cycle regulator at the G1/S transition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]