

Product datasheet for RC218427

Decorin (DCN) (NM_133506) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Decorin (DCN) (NM_133506) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: Decorin
Synonyms: CSCD; DSPG2; PG40; PGII; PGS2; SLRR1B
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC218427 representing NM_133506
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGGCCACTATCATCCTCCTTCTGCTTGCACAAGTTTCTGGGCTGGACCGTTTCAACAGAGAGGCT
 TATTTGACTTTATGCTAGAAGATGAGGCTTCTGGGATAGGCCAGAAAGTTCTGATGACCGGACTTCGA
 GCCCTCCCTAGGCCAGTGTCCTTCCGCTGTCAATGCCATCTTCGAGTGGTCCAGTGTCTGATTTG
 GGTCTGGACAAAGTGCCAAAGGATCTTCCCTGACACAACCTCTGCTAGACCTGCAAAACAACAAATAT
 CCGAAATCAAAGATGGAGACTTTAAGAACCTGAAGAACCTTACGTTGTCTACCTCATAACAACAATAT
 CTCTGTAGTTGGATCAAGTGACTTCTGCCACCTGGACACAACACCAAAAGGCTTCTTATTCGGGTGTG
 AGTCTTTTCAGCAACCCGGTCCAGTACTGGGAGATACAGCCATCCACCTCAGATGTGTCTACGTGCGCT
 CTGCCATTCAACTCGAAACTATAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_133506

ORF Size: 516 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_133506.3](#)


RefSeq Size: 1336 bp
RefSeq ORF: 519 bp
Locus ID: 1634
UniProt ID: [P07585](#)
Cytogenetics: 12q21.33
Domains: LRRNT, LRR
Protein Families: Druggable Genome, Secreted Protein
Protein Pathways: TGF-beta signaling pathway
MW: 17.4 kDa

Gene Summary: This gene encodes a member of the small leucine-rich proteoglycan family of proteins. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature protein. This protein plays a role in collagen fibril assembly. Binding of this protein to multiple cell surface receptors mediates its role in tumor suppression, including a stimulatory effect on autophagy and inflammation and an inhibitory effect on angiogenesis and tumorigenesis. This gene and the related gene biglycan are thought to be the result of a gene duplication. Mutations in this gene are associated with congenital stromal corneal dystrophy in human patients. [provided by RefSeq, Nov 2015]

Product images:



Circular map for RC218427



Coomassie blue staining of purified DCN protein (Cat# [TP318427]). The protein was produced from HEK293T cells transfected with DCN cDNA clone (Cat# RC218427) using MegaTran 2.0 (Cat# [TT210002]).