

Product datasheet for **RC218451L3V**

MAP2 (NM_031847) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MAP2 (NM_031847) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAP2
Synonyms:	MAP-2; MAP2A; MAP2B; MAP2C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031847
ORF Size:	1506 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218451).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031847.1
RefSeq Size:	2559 bp
RefSeq ORF:	1509 bp
Locus ID:	4133
UniProt ID:	P11137
Cytogenetics:	2q34
Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS
MW:	52.8 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described. [provided by RefSeq, Jan 2010]