

## Product datasheet for **RC218471L4V**

### **Munc18 1 (STXBP1) (NM\_003165) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	Munc18 1 (STXBP1) (NM_003165) Human Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	Munc18 1
<b>Synonyms:</b>	DEE4; MUNC18-1; N-Sec1; NSEC1; P67; RBSEC1; unc-18A; UNC18; unc18-1
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-mGFP-P2A-Puro (PS100093)
<b>Tag:</b>	mGFP
<b>ACCN:</b>	NM_003165
<b>ORF Size:</b>	1809 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(RC218471).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_003165.1</a>
<b>RefSeq Size:</b>	3899 bp
<b>RefSeq ORF:</b>	1812 bp
<b>Locus ID:</b>	6812
<b>UniProt ID:</b>	<a href="#">P61764</a>
<b>Cytogenetics:</b>	9q34.11
<b>Domains:</b>	Sec1
<b>MW:</b>	68.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a syntaxin-binding protein. The encoded protein appears to play a role in release of neurotransmitters via regulation of syntaxin, a transmembrane attachment protein receptor. Mutations in this gene have been associated with infantile epileptic encephalopathy-4. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2010]