

Product datasheet for **RC218473L3V**

TEM8 (ANTXR1) (NM_053034) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TEM8 (ANTXR1) (NM_053034) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TEM8
Synonyms:	ATR; GAPO; TEM8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053034
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218473).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053034.1
RefSeq Size:	1454 bp
RefSeq ORF:	1107 bp
Locus ID:	84168
UniProt ID:	Q9H6X2
Cytogenetics:	2p13.3
Domains:	VWA
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 38.3 kDa

Gene Summary: This gene encodes a type I transmembrane protein and is a tumor-specific endothelial marker that has been implicated in colorectal cancer. The encoded protein has been shown to also be a docking protein or receptor for Bacillus anthracis toxin, the causative agent of the disease, anthrax. The binding of the protective antigen (PA) component, of the tripartite anthrax toxin, to this receptor protein mediates delivery of toxin components to the cytosol of cells. Once inside the cell, the other two components of anthrax toxin, edema factor (EF) and lethal factor (LF) disrupt normal cellular processes. Three alternatively spliced variants that encode different protein isoforms have been described. [provided by RefSeq, Oct 2008]