

## Product datasheet for RC218492

### KRT13 (NM\_002274) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	KRT13 (NM_002274) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KRT13
Synonyms:	CK13; K13; WSN2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC218492 representing NM_002274 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGCCTCCGCCTGCAGAGCTCCTCTGCCAGCTATGGAGGTGGTTTCGGGGTGGCTCTTGCCAGCTGG  
GAGGAGGCCGTGGTGTCTCTACCTGTTCAACTCGGTTTGTGTCTGGGGATCAGCTGGGGCTATGGAGG  
CGGCGTGAGCTGTGGTTTGGTGGAGGGCTGGTGTAGTGGCTTTGGAGGTGGCTATGGAGGTGGCCTTGA  
GGTGGCTATGGAGGTGGCCTGGAGGTGGCTTTGGTGGGGTTTTGCTGGTGGCTTTGTTGACTTTGGT  
CTTGATGGCGCCTCCTCACTGGCAATGAGAAGATCACCATGCAGAACCTCAACGACCGCCTGGCTTC  
CTACCTGGAGAAGGTGCGCGCCCTGGAGGAGGCCAACGCTGACCTGGAGGTGAAGATCCGTGACTGGCAC  
CTGAAGCAGAGCCAGCTAGCCCTGAGCGGACTACAGCCCTACTACAAGACCATTGAAGAGCTCCGGG  
ACAAGATCCTGACCGCCACCATTTGAAAACAACCGGTCATCCTGGAGATTGACAATGCCAGGCTGGCTGC  
GGACGACTTCAGGCTCAAGTATGAGAAATGAGCTGGCCCTGCGCCAGAGCGTGGAGGCCGACATCAACGGC  
CTGCGCCGGTGTGGATGAGCTCACTCTGTCTAAGACTGACCTGGAGATGCAGATCGAGAGCCTGAATG  
AAGAGCTAGCCTACATGAAGAAGAACCATGAAGAGGAGATGAAGGAATTTAGCAACCAGGTGGTCCGCCA  
GGTCAACGTGGAGATGGATGCCACCCAGGCATTGACCTGACCCGCTGCTGGCAGAGATGAGGGAGCAG  
TACGAGCCATGGCAGAGAGGAACCGCGGATGCTGAGGAATGGTTCCACGCCAAGAGTGCAGAGCTGA  
ACAAGGAGGTGTCTACCAACACTGCCATGATTCAGACCAGCAAGACAGAGATCACGGAGCTCAGGCGCAC  
GCTCCAAGGCCTGGAGATTGAGCTGCAGTCCCAGCTGAGCATGAAAGCGGGGCTGGAGAACCGGTGGCA  
GAGACGGAGTGCCGCTATGCCCTGCAGCTGCAGCAGATCCAGGGACTCATCAGCAGCATCGAGGCCAGC  
TGAGCGAGCTCCGAGTGAGATGGAGTGCAGAACCAAGAGTACAAGATGCTGCTGGACATCAAGACAGC  
TCTGGAGCAGGAGATCGCCACCTACCGCAGCCTGCTCGAGGGCCAGGACGCCAAGAAGCGTCAGCCCCG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



**Protein Sequence:** >RC218492 representing NM\_002274  
 Red=Cloning site Green=Tags(s)



MSLR LQSSSASYGGGFGGGSCQLGGGRGVSTCSTRFVSGGSAGGYGGVSCGFGGGAGSGFGGGYGGGLG  
 GGYGGGLGGGFGGFAGGFVDFGACDGLLTGNEKITMQNLNDR LASYLEKVRALEEANADLEVKIRDWH  
 LKQSPASPERDYSPYYKTEELRDKILTATIENNRVILEIDNARLAADDFRLKYENELALRQSVEADING  
 LRRVLDLTL SKTDLEMQIESLNEELAYMKKNHEEEMKEFSNQVVGQVNVEMDATPGIDL TRVLAEMREQ  
 YEAMAERNRRDAEEWFHAKSAELNKEVSTNTAMIQT SKTEITELRRTLQGLEIELQSQLSMKAGLENTVA  
 ETECRYALQLQQIQLISSIEAQLSELRSEMECQNQEYKMLLDIKTRLEQEIATYRSLLLEGQDAKKRQPP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8058\\_f04.zip](https://cdn.origene.com/chromatograms/mk8058_f04.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_002274

**ORF Size:** 1260 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_002274.3](#)

**RefSeq Size:** 1693 bp

**RefSeq ORF:** 1263 bp

**Locus ID:** 3860

**UniProt ID:** [P13646](#)


**Cytogenetics:** 17q21.2

**Domains:** filament

**MW:** 45.7 kDa

**Gene Summary:** The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This type I cytokeratin is paired with keratin 4 and expressed in the suprabasal layers of non-cornified stratified epithelia. Mutations in this gene and keratin 4 have been associated with the autosomal dominant disorder White Sponge Nevus. The type I cytokeratins are clustered in a region of chromosome 17q21.2. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been described. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for RC218492