

## Product datasheet for **RC218498L4V**

### Syntabulin (SYBU) (NM\_001099756) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Syntabulin (SYBU) (NM_001099756) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Syntabulin
Synonyms:	GOLSYN; OCSYN; SNPHL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001099756
ORF Size:	1980 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218498).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001099756.1</a> , <a href="#">NP_001093226.1</a>
RefSeq Size:	2969 bp
RefSeq ORF:	1983 bp
Locus ID:	55638
UniProt ID:	<a href="#">Q9NX95</a>
Cytogenetics:	8q23.2
MW:	71.9 kDa



[View online »](#)

**Gene Summary:**

Syntabulin/GOLSYN is part of a kinesin motor-adaptor complex that is critical for the anterograde axonal transport of active zone components and contributes to activity-dependent presynaptic assembly during neuronal development (Cai et al., 2007 [PubMed 17611281]).[supplied by OMIM, Mar 2008]