

Product datasheet for **RC218522**

C8orf58 (NM_001013842) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C8orf58 (NM_001013842) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C8orf58
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC218522 representing NM_001013842 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGGGCCGGCGCGCCTTCGCCGTGGACGGCCGGGATGGGGCTGGCGAGGGCCTGGCACGGGGCT
GCATAGTGCTGGAGTCACCGACCTACAGACGGATCCCCGACGCTGCCACGGGTGCTCATCTGGGA
GAGAGGTGACAAGTTCAGAGGTGTCGGCAGGGAGGCACTCTTTCTCAAAGTGGCTCCCGGGACTCAGGA
GTGGAGATGGCAGTTGGGGACAGCCCCCTGGCCGCCTTACCGGGCCTTTCTCAGGACTCCCTGGACTTTG
AATCCTCAGGGAGTTCTGAGCCCCCGCCAGGTAGGCCGACTCCTGGCCAGCCAGAAGCTGGGGGAGGT
GTTGGAGCGGTCCCGCCGGCTCCCAACAGCTCCACAGCTTGTCAGGACAACACCGCTCCCTGCGGCTG
GCAAGCAAGCCTGAGCGTGAAGTGCCCTTGAGCAGGGCAACAGGAGTCCATGGAGGCAGACACAGACC
TAGAGGCAGGCCTGGAAGAAGAGGCGGTGGGGGCCTGGGGCCTGGAGCCTGGGCTGCCTCCCCGGGCA
GGGTCTTCGCTATCTGGAACACCTGTGCCTGGTGTGCTGGAGCAGATGGCAAGGCTCCAGCAGCTCTACCTG
CAGCTGCGGATCCAGAGGCCCCAGGGGATCCCGGCGAGGAGGAGTCGACCCGAGCCCCCTTACCGTCCC
CGTTACACACCCAGGCAATCGGGGCAGGGGCCATGGGAGCTGCTAAGCCAGACAGAGCACACAGGAGC
AAAGGCTGCTTACCCCCAAAGGTGGAGGTGCCAGTGCCAACCTCCAGGCTGCCAGAAACCCAGTG
GAGCCAACGTACCACTTGCCATCTCCAGGGACACAAGCGGGATATCTCCACTGGGACAAGGTCAAGG
TCTGTCTAACCGGATCTGCCGAGAAGCCACCACCCTGAGCCCCCTGCCCTCCTGATGGCTCTGA
CCCCAGGATCGAGTCCAGGGACCTCCCTGAAAGGCCTCAGTGCCGCCCCACCGGAAGACCTTTATGCCA
TCATTAGTGGTTAAGAAGCAACGAGCAAAAAACCTTTCTGTAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC218522 representing NM_001013842
 Red=Cloning site Green=Tags(s)


MMGRRRAFAVDGRDGAGEGLARGCIVPGVTSTYRRIPDAAHGCSWGERGDKFRGVGREALFLKLASRDG
 VEMAVGDSPLAALPGLSQDSLDFESSGSSEPPAQVGRLLASQKLGEVLERSRRLPTAPTSLSGQHRSLRL
 ASKPEREVPLGAGQQESMEADTDLEAGLEEEAVGGLGPGAWACLPQGGLRYLEHLCLVLEQMARLQQLYL
 QLRIQRPPGDPGEEESTRAPLPSPLHTPGNRGQGPWELLSQTEHTGAKAASPPKVEVPSANPPRLPETPV
 EPTYHLPSQGHKRDISHWDKVKVLLNRRICRRSHHHPEPPAPDGDSPRIESRDLPERPQCRPHRKTFFMP
 SLVVKQRAKNLSVG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8021_f12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001013842

ORF Size: 1095 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001013842.3

RefSeq Size: 2075 bp

RefSeq ORF: 1098 bp


Locus ID: 541565

UniProt ID: Q8NAV2


Cytogenetics: 8p21.3

MW: 39.5 kDa


Product images:



Circular map for RC218522



Western blot validation of overexpression lysate (Cat# [LY423043]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC218522 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified C8orf58 protein (Cat# [TP318522]). The protein was produced from HEK293T cells transfected with C8orf58 cDNA clone (Cat# RC218522) using MegaTran 2.0 (Cat# [TT210002]).