

## Product datasheet for **RC218530L4V**

### **T (NM\_003181) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	T (NM_003181) Human Tagged ORF Clone Lentiviral Particle
Symbol:	T
Synonyms:	SAVA; T; TFT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003181
ORF Size:	1308 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218530).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003181.2</a>
RefSeq Size:	2518 bp
RefSeq ORF:	1308 bp
Locus ID:	6862
UniProt ID:	<a href="#">O15178</a>
Cytogenetics:	6q27
Protein Families:	ES Cell Differentiation/IPS, Transcription Factors
MW:	47.3 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is an embryonic nuclear transcription factor that binds to a specific DNA element, the palindromic T-site. It binds through a region in its N-terminus, called the T-box, and effects transcription of genes required for mesoderm formation and differentiation. The protein is localized to notochord-derived cells. Variation in this gene was associated with susceptibility to neural tube defects and chordoma. A mutation in this gene was found in a family with sacral agenesis with vertebral anomalies. [provided by RefSeq, Sep 2018]