

Product datasheet for **RC218541L3V**

DUOXA2 (NM_207581) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | DUOXA2 (NM_207581) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | DUOXA2 |
| Synonyms: | SIMNIPHOM; TDH5 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_207581 |
| ORF Size: | 960 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC218541). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_207581.2 |
| RefSeq Size: | 1440 bp |
| RefSeq ORF: | 963 bp |
| Locus ID: | 405753 |
| UniProt ID: | Q1HG44 |
| Cytogenetics: | 15q21.1 |
| Protein Families: | Transmembrane |
| MW: | 34.6 kDa |



[View online »](#)

Gene Summary:

This gene encodes an endoplasmic reticulum protein that is necessary for proper cellular localization and maturation of functional dual oxidase 2. Mutations in this gene have been associated with thyroid dyshormonogenesis 5.[provided by RefSeq, Feb 2010]