

Product datasheet for **RC218645A1V**

Human TCEA1 (NM_201437) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TCEA1 (NM_201437) AAV Particle
Tag:	Myc-DDK
Symbol:	TCEA1
Synonyms:	GTF2S; SII; TCEA; TF2S; TFII5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC218645 representing NM_201437 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGGACGAAGTGGTCCGCTTTGCCAAGAAGATGGACAAGATGGTGCAGAAGAAGACGCGTCCACAA
GAATCGGAATGTCAGTTAATGCTATTCGCAAGCAGAGTACAGATGAGGAAGTTACATCTTTGGCAAAGTC
TCTCATCAAATCCTGGAAAAATTATTAGATGGGCCATCAACTGAGAAAGACCTTGACGAAAAGAAGAAA
GAACCTGCAATTACATCGCAGAACAGCCCTGAGGCAAGAGAAGAAAGTACTTCCAGCGGCAATGTAAGCA
ACAGAAAGGATGAGACAAATGCTCGAGATACTTATGTTTCATCCTTTCTCGGGCACCAAGCACTTCTGA
TTCTGTGCGGTTGAAGTGTAGGGAGATGCTTGCTGCAGCTCTTCGAACAGGGGATGACTACATTGCAATT
GGAGCTGATGAGGAAGATTAGGATCTCAAATTGAAGAAGCTATATATCAAGAAATAAGGAATACAGACA
TGAAATACAAAAATAGAGTACGAAGTAGGATATCAAATCTTAAAGATGCAAAAAATCCAAATTTAAGGAA
AAATGTCCTCTGTGGGAATATTCTCTGACTTATTTGCTAGAATGACAGCAGAGGAAATGGCTAGTGAT
GAGCTGAAAGAGATGCGGAAAACTTGACCAAGAAGCCATCAGAGAGCATCAGATGGCCAAGACTGGTG
GGACCCAGACTGACTTGTTACATGTGGCAAATGTAAAAAGAAGATTGCACTTACACACAGGTACAAAC
CCGTAGTGCTGATGAACCAATGACAACATTTGTTGTCTGTAATGAATGTGGAATCGATGGAAGTTCTGT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC218645 representing NM_201437 Red=Cloning site Green=Tags(s) MEDEVVRFakkMDKmvQkKNAstrIGMSVNAIRKQSTDEEVtSLAKSLIKSWKKLLDGPSTeKDLDEKKK EPAITsQNSPEAREESTSSGNVSNRKDEtNARDTYVSSFPRAPSTSDSVRLKCREMLAAALRTGDDYIAI GADEEELGSQIEEAIYQEIRNTDMKYKNRVRSRISNLKDAKNPNLRKNVLCGNIPPDLFARMTAEEMASD ELKEMRKNLTKEAIREHQMAKTGGTQTDLFTCGKCKKKNCTYTQVQTRSADEPMTTFVVCNECGNRWKFC SGPTrTRRLEqKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_201437
ORF Size:	840 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_201437.1
RefSeq Size:	2721 bp
RefSeq ORF:	843 bp
Locus ID:	6917
UniProt ID:	P23193
Cytogenetics:	8q11.23
MW:	31.5 kDa