

Product datasheet for RC218671

OR9K2 (NM_001005243) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR9K2 (NM_001005243) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR9K2
Synonyms:	OR12-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC218671 representing NM_001005243 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTAGGATCCAAACCAAGAGTTCATTTGTATATTTGCCCTGTGCCTCTCAACAGTTTCTACCATGG
GTGACAGGGGAACAAGCAATCACTCAGAAATGACTGACTTCATTCTTGACGGCTTCAGGGTACGCCCAGA
GCTCCACATTCTCCTCTCCTGCTATTTTGTGGTTATGCCATGATCCTTCTAGGGAATGTTGGGATG
ATGACCATTATTGACTGATCCTCGGCTGAACACACCAATGTATTTTCTAGGCAATCTCCTCTCA
TTGATCTTTTCTATTCATCTGTTATTGAACCAAGGCTATGATCAACTTCTGGTCTGAAAACAAGTCTAT
CTCCTTTGACGGCTGTGTGGCCAGCTTTTCTTTGCCCTCCTCATTGTGACTGAGGGATTTCTCCTG
GCGGCCATGGCTTATGACCGCTTATTGCCATCTGCAACCCTCTGCTCTACTCTGTTCAAATGTCCACAC
GTCTGTGTAAGTCAAGTGGTGGCTGGTTCCTATTTTGTGGCTGCATTAGCTCAGTTATTCAGACTAGCAT
GACATTTACTTTATCTTTTGGCCTTCTCGGGCTGTTGACCACCTTTACTGTGATTCTCGCCCACTTCAG
AGACTGTCTTGTCTGATCTCTTATCCATAGAATGATATCTTTTCTTATCATGTATTATTATCTTGC
CTACTATCATAGTCATTATAGTATCTTACATGTATTTGTGTCCACAGTTCTAAAGATACATTCTACTGA
GGGACATAAGAAGGCCTTCTCCACCTGCAGCTCTCACCTGGGAGTTGTGAGTGTGCTGTATGGTGTCTGC
TTTTTTATGTATCTCACTCCTGACAGATTTCTGAGCTGAGTAAAGTGGCATCCTTATGTTACTCCCTAG
TCACTCCATGTTGAATCCTTTGATTTACTCTCTGAGGAACAAAGATGTCCAAGAGGCTCTAAAAAATT
TCTAGAGAAGAAAAATATTATTCTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC218671 representing NM_001005243
Red=Cloning site Green=Tags(s)

MLGSKPRVHLYILPCASQQVSTMGDRGTSNHSEMTDFILAGFRVRPELHILLFLLFLFVYAMILLGNVGM
 MTIIMTDPRLNTPMYFFLGNLSFIDLFYSSVIEPKAMINFWSENKISIFAGCVAQLFLFALLIVTEGFL
 AAMAYDRFIAICNPLLYSQMSTRCTQLVAGSYFCGCISSVIQTSMTFTLSFCASRAVDHFYCDSPRLQ
 RLSCSDLFIHRMISFSLSCIIILPTIIIVIVSYMYIVSTVLKIHSTEGHKKAFSTCSSHLGVVSVLYGAV
 FFMYLTPDRFPELSKVASLCYSLVTPMLNPLIYSLRNKDVQEALKKFLEKKNIIIL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8003_e12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001005243

ORF Size: 1005 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005243.1](#), [NP_001005243.1](#)

RefSeq Size: 1008 bp

RefSeq ORF: 942 bp

Locus ID: 441639

UniProt ID: [Q8NGE7](#)

Cytogenetics: 12q13.2


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 37.6 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC218671