

## Product datasheet for **RC218676L2V**

### **AFAP1L1 (NM\_152406) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	AFAP1L1 (NM_152406) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AFAP1L1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_152406
ORF Size:	2304 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218676).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152406.1</a>
RefSeq Size:	4223 bp
RefSeq ORF:	2307 bp
Locus ID:	134265
UniProt ID:	<a href="#">Q8TED9</a>
Cytogenetics:	5q32
MW:	86.4 kDa
Gene Summary:	May be involved in podosome and invadosome formation.[UniProtKB/Swiss-Prot Function]


[View online »](#)