

Product datasheet for **RC218678A1V**

Human INMT (NM_006774) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human INMT (NM_006774) AAV Particle
Tag:	Myc-DDK
Symbol:	INMT
Synonyms:	TEMT
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC218678 representing NM_006774 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAGGGTGGCTTCACTGGGGGTGATGAGTACCAGAAGCACTTCCTGCCAGGGACTACTTGGCTACTT
ACTACAGCTTCAATGGCAGCCCTCACCCGAGGCCGAGATGCTGAAGTTAACTTGAATGTCTCCACAA
GACCTTCGGCCCTGGAGGCCTCCAAGGGGACACGCTGATTGACATTGGCTCAGGTCTACCATCTACCAA
GTTCTTGCTGCCTGTGATTCCTTCCAAGACATCACTCTCTCCGACTTTACCGACCGCAACCGGGAGGAGC
TGGAAAAGTGGCTGAAGAAGGAGCCGGGGGCTATGACTGGACCCAGCGGTGAAATTCGCCTGTGAGCT
GGAAGGAAACAGCGGCCGATGGGAGGAGAAGGAGGAGAAGCTGCGGGCAGCGGTGAAGCGGGTGTCTCAAG
TGCGATGTCCACCTGGGCAACCCGCTGGCCCGGCTGTGTTGCCTCTCGCCGACTGTGTGCTCACCTGC
TGGCCATGGAGTGTGCCTGCTGTAGCCTTGATGCCTACCGCGCTGCCCTGTGCAACCTTGCCCTACTGCT
CAAGCCGGGTGGCCACCTGGTGACCACTGTACGCTTCGGCTCCCGTCTACGTGGTGGGAAGCGTGAA
TTTCTGCTGGCCCTGGAGAAAGAGGAGGTGGAGCAGGCTGCTCGGATGCTGGCTTTGACATTGAAC
AGCTCTACACAGTCCCAGAGCTACTGTGTCACCAATGCTGCCAACAATGGGGTCTGCTGCATTGTGGC
TCGCAAGAAGCTGGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC218678 representing NM_006774 Red=Cloning site Green=Tags(s) MKGGFTGGDEYQKHFLPRDYLATYYSFNGSPSPEAEMLKFNLECLHKTFGPGGLQGDTLIDIGSGPTIYQ VLAACDSFQDITLSDFTDNRNREELEKWLKKEPGAYDWTAPVKFACELEGNSEGRWEEKKELRAAVKRVLK CDVHLGNPLAPAVLPLADCVLTLLAMECACCSLDAYRAALCNLASLLKPGGHLVTTVTLRPSYVVGKRE FSCVALEKEEVEQAVLDAGFDIEQLLHSPQSYSVTNAANNGVCCIVARKKPGP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006774
ORF Size:	789 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006774.4</u>
RefSeq Size:	2639 bp
RefSeq ORF:	792 bp
Locus ID:	11185
UniProt ID:	<u>O95050</u>
Cytogenetics:	7p14.3
MW:	28.7 kDa