

## Product datasheet for **RC218682L1V**

### TRAF3 (NM\_003300) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | TRAF3 (NM_003300) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | TRAF3  |
| Synonyms:                 | CAP-1; CAP1; CD40bp; CRAF1; IIAE5; LAP1; RNF118  |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)  |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_003300  |
| ORF Size:                 | 1704 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC218682).   |
| OTI Disclaimer:           | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_003300.2</a>  |
| RefSeq Size:              | 7654 bp  |
| RefSeq ORF:               | 1707 bp  |



[View online »](#)

|                          |   |
|--------------------------|---|
| <b>Locus ID:</b>         | 7187  |
| <b>UniProt ID:</b>       | <a href="#">Q13114</a>  |
| <b>Cytogenetics:</b>     | 14q32.32  |
| <b>Domains:</b>          | zf-TRAF, RING, MATH   |
| <b>Protein Families:</b> | Druggable Genome  |
| <b>Protein Pathways:</b> | Pathways in cancer, RIG-I-like receptor signaling pathway, Small cell lung cancer, Toll-like receptor signaling pathway   |
| <b>MW:</b>               | 64.5 kDa  |
| <b>Gene Summary:</b>     | <p>The protein encoded by this gene is a member of the TNF receptor associated factor (TRAF) protein family. TRAF proteins associate with, and mediate the signal transduction from, members of the TNF receptor (TNFR) superfamily. This protein participates in the signal transduction of CD40, a TNFR family member important for the activation of the immune response. This protein is found to be a critical component of the lymphotoxin-beta receptor (LTbetaR) signaling complex, which induces NF-kappaB activation and cell death initiated by LTbeta ligation. Epstein-Barr virus encoded latent infection membrane protein-1 (LMP1) can interact with this and several other members of the TRAF family, which may be essential for the oncogenic effects of LMP1. The protein also plays a role in the regulation of antiviral response. Mutations in this are associated with Encephalopathy, acute, infection-induced, herpes-specific 5. [provided by RefSeq, Jul 2020]</p> |