

Product datasheet for RC218713

H2AZ2 (NM_201517) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2AZ2 (NM_201517) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: H2AZ2
Synonyms: H2A.Z-2; H2AFV; H2AV
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC218713 representing NM_201517
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGCTGGAGGCAAAGCTGGAAGGACAGTGGGAAGGCCAAGGCTAAGGCAGTATCTCGCTCACAGAGAG
 CTGGGCTACAGGTGCTGGAGCTGGCAGGTAATGCTTCTAAGGATCTCAAAGTAAAGCGTATCACTCCGCG
 TCACTTGACGCTTGCAATCCGTGGTGAAGAGTTGGATTCTCTATCAAGGCTACCATAGCTGGGGT
 GGTGTGATCCCTCACATCCACAAATCTCTGATTGAAAGAAGGGACAGCAGAAAAGTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC218713 representing NM_201517
 Red=Cloning site Green=Tags(s)

MAGGKAGKDSGKAKAKAVSRSQRAGLQVLELAGNASKDLKVKRITPRHLQLAIRGDEELDSLKATIAGG
 GVPHIHKSLIGKKGQKTA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_201517

ORF Size: 270 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_201517.2</u> , <u>NP_958925.1</u>
RefSeq Size:	1315 bp
RefSeq ORF:	273 bp
Locus ID:	94239
UniProt ID:	<u>Q71UI9</u>
Cytogenetics:	7p13
Protein Families:	Druggable Genome
Protein Pathways:	Systemic lupus erythematosus
MW:	9.4 kDa
Gene Summary:	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent histone that is a member of the histone H2A family. Several transcript variants encoding different isoforms, have been identified for this gene. [provided by RefSeq, Oct 2015]