

Product datasheet for **RC218770L2V**

STK19 (NM_004197) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	STK19 (NM_004197) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STK19
Synonyms:	D6S60; D6S60E; G11; HLA-RP1; RP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004197
ORF Size:	1092 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218770).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004197.1 , NP_004188.1
RefSeq Size:	1620 bp
RefSeq ORF:	765 bp
Locus ID:	8859
UniProt ID:	P49842
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome, Protein Kinase
MW:	40.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a serine/threonine kinase which localizes predominantly to the nucleus. Its specific function is unknown; it is possible that phosphorylation of this protein is involved in transcriptional regulation. This gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6 and expresses two transcript variants. [provided by RefSeq, Jul 2008]