

Product datasheet for **RC218787L4V**

Cell adhesion molecule 2 (CADM2) (NM_153184) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cell adhesion molecule 2 (CADM2) (NM_153184) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cell adhesion molecule 2
Synonyms:	IGSF4D; Necl-3; NECL3; SynCAM 2; synCAM2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153184
ORF Size:	1311 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218787).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153184.2 , NP_694854.2
RefSeq Size:	3347 bp
RefSeq ORF:	1314 bp
Locus ID:	253559
UniProt ID:	Q8N3J6
Cytogenetics:	3p12.1
Domains:	ig, IGc2, IG
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 47.4 kDa

Gene Summary: This gene encodes a member of the synaptic cell adhesion molecule 1 (SynCAM) family which belongs to the immunoglobulin (Ig) superfamily. The encoded protein has three Ig-like domains and a cytosolic protein 4.1 binding site near the C-terminus. Proteins belonging to the protein 4.1 family crosslink spectrin and interact with other cytoskeletal proteins. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]