

Product datasheet for RC218813L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD47 (NM_001777) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: CD47 (NM 001777) Human Tagged ORF Clone Lentiviral Particle

Symbol: CD47

Synonyms: IAP; MER6; OA3

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK
ACCN: NM 001777

ORF Size: 969 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC218813).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 001777.3

 RefSeq Size:
 5346 bp

 RefSeq ORF:
 972 bp

 Locus ID:
 961

 UniProt ID:
 Q08722

Cytogenetics: 3q13.12

Domains: CD47

Protein Families: Druggable Genome, Transmembrane





CD47 (NM_001777) Human Tagged ORF Clone Lentiviral Particle - RC218813L3V

Protein Pathways: ECM-receptor interaction

MW: 35.21 kDa

Gene Summary: This gene encodes a membrane protein, which is involved in the increase in intracellular

calcium concentration that occurs upon cell adhesion to extracellular matrix. The encoded protein is also a receptor for the C-terminal cell binding domain of thrombospondin, and it may play a role in membrane transport and signal transduction. This gene has broad tissue distribution, and is reduced in expression on Rh erythrocytes. Alternatively spliced transcript

variants have been found for this gene. [provided by RefSeq, Jul 2010]