

Product datasheet for **RC218846L4V**

BCKDHB (NM_000056) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BCKDHB (NM_000056) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BCKDHB
Synonyms:	BCKDE1B; BCKDH E1-beta; E1B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000056
ORF Size:	1176 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218846).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000056.2
RefSeq Size:	1552 bp
RefSeq ORF:	1179 bp
Locus ID:	594
UniProt ID:	P21953
Cytogenetics:	6q14.1
Domains:	transket_pyr, transketolase_C
Protein Pathways:	Metabolic pathways, Valine, leucine and isoleucine degradation



[View online »](#)

MW: 43.12 kDa

Gene Summary: This gene encodes the E1 beta subunit of branched-chain keto acid dehydrogenase, which is a multienzyme complex associated with the inner membrane of mitochondria. This enzyme complex functions in the catabolism of branched-chain amino acids. Mutations in this gene have been associated with maple syrup urine disease (MSUD), type 1B, a disease characterized by a maple syrup odor to the urine in addition to mental and physical retardation and feeding problems. Alternative splicing at this locus results in multiple transcript variants. [provided by RefSeq, Jan 2016]