

Product datasheet for **RC218849L1V**

EGFL7 (NM_016215) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EGFL7 (NM_016215) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EGFL7
Synonyms:	NEU1; VE-STATIN; ZNEU1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016215
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218849).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016215.3
RefSeq Size:	1569 bp
RefSeq ORF:	822 bp
Locus ID:	51162
UniProt ID:	Q9UHF1
Cytogenetics:	9q34.3
Domains:	EGF_CA, EGF, EGF
Protein Families:	Secreted Protein



[View online »](#)

MW: 29.4 kDa

Gene Summary: This gene encodes a secreted endothelial cell protein that contains two epidermal growth factor-like domains. The encoded protein may play a role in regulating vasculogenesis. This protein may be involved in the growth and proliferation of tumor cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]