

Product datasheet for **RC218860L3V**

CDRT1 (NM_006382) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CDRT1 (NM_006382) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDRT1
Synonyms:	C17ORF1; C17ORF1A; C17ORF1; FBXW10B; FBXW10P1; HREP; SM25H2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006382
ORF Size:	2262 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218860).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006382.3
RefSeq Size:	2780 bp
RefSeq ORF:	2259 bp
Locus ID:	374286
UniProt ID:	O95170
Cytogenetics:	17p12
Domains:	WD40
MW:	85.6 kDa



[View online »](#)

Gene Summary:

Members of the F-box protein family, such as FBXW10, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603034), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]