

Product datasheet for **RC218864L1V**

LRP5 (NM_002335) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LRP5 (NM_002335) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LRP5
Synonyms:	BMND1; EVR1; EVR4; HBM; LR3; LRP-5; LRP-7; LRP7; OPPG; OPS; OPTA1; PCLD4; VBCH2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002335
ORF Size:	4845 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218864).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002335.1
RefSeq Size:	5100 bp
RefSeq ORF:	4848 bp
Locus ID:	4041
UniProt ID:	O75197
Cytogenetics:	11q13.2
Domains:	Idl_recept_b, Idl_recept_a, EGF, EGF
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: Wnt signaling pathway

MW: 179 kDa

Gene Summary: This gene encodes a transmembrane low-density lipoprotein receptor that binds and internalizes ligands in the process of receptor-mediated endocytosis. This protein also acts as a co-receptor with Frizzled protein family members for transducing signals by Wnt proteins and was originally cloned on the basis of its association with type 1 diabetes mellitus in humans. This protein plays a key role in skeletal homeostasis and many bone density related diseases are caused by mutations in this gene. Mutations in this gene also cause familial exudative vitreoretinopathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014]