

Product datasheet for **RC218886L3V**

INPP4A (NM_001566) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	INPP4A (NM_001566) Human Tagged ORF Clone Lentiviral Particle
Symbol:	INPP4A
Synonyms:	INPP4; TVAS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001566
ORF Size:	2862 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218886).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001566.1 , NP_001557.1
RefSeq Size:	3224 bp
RefSeq ORF:	2865 bp
Locus ID:	3631
UniProt ID:	Q96PE3
Cytogenetics:	2q11.2
Protein Families:	Transmembrane
Protein Pathways:	Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system



[View online »](#)

MW: 107.1 kDa

Gene Summary: This gene encodes an Mg⁺⁺ independent enzyme that hydrolyzes the 4-position phosphate from the inositol ring of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate, and inositol 3,4-bisphosphate. Multiple transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Aug 2008]