

Product datasheet for **RC218896L1V**

VGLL3 (NM_016206) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VGLL3 (NM_016206) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VGLL3
Synonyms:	VGL-3; VGL3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016206
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218896).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016206.2
RefSeq Size:	10396 bp
RefSeq ORF:	981 bp
Locus ID:	389136
UniProt ID:	A8MV65
Cytogenetics:	3p12.1
Protein Families:	Transcription Factors
MW:	36 kDa



[View online »](#)

Gene Summary: May act as a specific coactivator for the mammalian TEFs.[UniProtKB/Swiss-Prot Function]