

Product datasheet for RC218907L3V

HERC6 (NM_017912) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HERC6 (NM_017912) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HERC6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017912
ORF Size:	3066 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218907).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017912.3
RefSeq Size:	3903 bp
RefSeq ORF:	3069 bp
Locus ID:	55008
UniProt ID:	Q8IVU3
Cytogenetics:	4q22.1
Domains:	HECT
Protein Families:	Druggable Genome
MW:	114.9 kDa


[View online »](#)

Gene Summary:

HERC6 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).[supplied by OMIM, Mar 2008]