

Product datasheet for **RC218912L4V**

FAM3C (NM_001040020) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | FAM3C (NM_001040020) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | FAM3C |
| Synonyms: | GS3786; ILEI |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001040020 |
| ORF Size: | 681 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC218912). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001040020.1 , NP_001035109.1 |
| RefSeq Size: | 2512 bp |
| RefSeq ORF: | 684 bp |
| Locus ID: | 10447 |
| UniProt ID: | Q92520 |
| Cytogenetics: | 7q31.31 |
| Protein Families: | Secreted Protein, Transmembrane |
| MW: | 24.7 kDa |



[View online »](#)

Gene Summary:

This gene is a member of the family with sequence similarity 3 (FAM3) family and encodes a secreted protein with a GG domain. A change in expression of this protein has been noted in pancreatic cancer-derived cells. [provided by RefSeq, Mar 2010]