

## Product datasheet for **RC218914A1V**

### Human S100 Calcium Binding Protein A13 (S100A13) (NM\_001024212) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human S100 Calcium Binding Protein A13 (S100A13) (NM\_001024212) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** S100 Calcium Binding Protein A13  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC218914 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCAGCAGAACCCTGACAGAGCTAGAGGAGTCCATTGAGACCGTGGTCACCACCTTCTTCACCTTTG  
CAAGGCAGGAGGGCCGGAAGGATAGCCTCAGCGTCAACGAGTTCAAAGAGCTGGTTACCCAGCAGTTGCC  
CCATCTGCTCAAGGATGTGGGCTCTCTTGATGAGAAGATGAAGAGCTGGATGTGAATCAGGACTCGGAG  
CTCAAGTTCAATGAGTACTGGAGATTGATTGGGGAGCTGGCCAAGGAAATCAGGAAGAAGAAAGACCTGA  
AGATCAGGAAGAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC218914 protein sequence  
Red=Cloning site Green=Tags(s)

MAAEPLTELEESIETVVTFFTFARQEGRKDSLVSNEFKELVTQQLPHLLKDVGSLDEKMKSLDVNQDSE  
LKFNEYWRLIGELAKEIRKKKDLKIRKK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001024212  
**ORF Size:** 294 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

**Stability:** AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

**RefSeq:** [NM\\_001024212.1](#), [NP\\_001019383.1](#)

**RefSeq Size:** 638 bp

**RefSeq ORF:** 297 bp

**Locus ID:** 6284

**UniProt ID:** [Q99584](#)

**Cytogenetics:** 1q21.3

**MW:** 11.5 kDa