

Product datasheet for **RC218920A1V**

Human Glucokinase (GCK) (NM_033508) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Glucokinase (GCK) (NM_033508) AAV Particle
Tag:	Myc-DDK
Symbol:	Glucokinase
Synonyms:	FGQTL3; GK; GLK; HHF3; HK4; HKIV; HXKP; LGLK; MODY2; PNDM1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

ORF Nucleotide Sequence:

>RC218920 representing NM_033508
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCAGACCAAGATCCCAACTCCACAACCCAACTCCAGGTAGAGCAGATCCTGGCAGAGTTCCAGC
 TGCAGGAGGAGGACCTGAAGAAGGTGATGAGACGGATGCAGAAGGAGATGGACCGCGCCTGAGGCTGGA
 GACCCATGAAGAGGCCAGTGTGAAGATGCTGCCACCTACGTGCGCTCCACCCAGAAGGCTCAGAAGTC
 GGGGACTTCTCTCCCTGGACCTGGGTGGCACTAACTTCAGGGTGTGCTGGTGAAGGTGGGAGAAGGTG
 AGGAGGGGAGTGGAGCGTGAAGACCAACACCAGATGTACTCCATCCCCGAGGACGCCATGACCGGCAC
 TGCTGAGATGCTCTCGACTACATCTCTGAGTGCATCTCCGACTTCTGGACAAGCATCAGATGAAACAC
 AAGAAGCTGCCCTGGGCTTACCTTCTCTTCTGTGAGGCACGAAGACATCGATAAGGGCATCTTTC
 TCAACTGGACCAAGGGCTTCAAGGCCTCAGGAGCAGAAGGGAACAATGTCGTGGGGCTTCTGCGAGACGC
 TATCAAACGGAGAGGGGACTTTGAAATGGATGTGGTGGCAATGGTGAATGACACGGTGGCCACGATGATC
 TCCTGCTACTACGAAGACCATCAGTGCAGGTCGGCATGATCGTGGGCACGGGCTGCAATGCCTGCTACA
 TGGAGGAGATGCAGAATGTGGAGCTGGTGGAGGGGACGAGGGCCGCATGTGCGTCAATACCGAGTGGGG
 CGCCTTCGGGGACTCCGGCGAGCTGGACGAGTTCTGCTGGAGTATGACCGCCTGGTGGACGAGAGCTCT
 GCAAACCCCGGTGAGCAGCTGTATGAGAAGCTCATAGTTGGCAAGTACATGGGCGAGCTGGTGGCGCTTG
 TGCTGCTCAGGCTCGTGACGAAAACCTGCTCTTCCACGGGGAGGCTCCGAGCAGCTGCGCACACGCGG
 AGCCTTCGAGACGCGCTTCGTGTGCGAGGTGGAGAGCGACACGGGCGACCGCAAGCAGATCTACAACATC
 CTGAGCACGCTGGGGCTGCGACCTCGACCACGACTGCGACATCGTGGCGCGCGCTGCGAGAGCGTGT
 CTACGCGCGCTGCGACATGTGCTCGGCGGGGCTGGCGGGCGTCATCAACCGCATGCGCGAGAGCCGAG
 CGAGGACGTAATGCGCATCACTGTGGCGTGGATGGCTCCGTGTACAAGCTGCACCCAGCTTCAAGGAG
 CGGTTCCATGCCAGCGTGCAGGCTGACGCCAGCTGCGAGATCACCTTCATCGAGTCGGAGGAGGGCA
 GTGGCCGGGGCGGGCCCTGGTCTCGCGGTGGCTGTGAAGAAGGCTGTATGCTGGGCCAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC218920 representing NM_033508
 Red=Cloning site Green=Tags(s)

MPRPRSQLPQNSQVEQILAEFQLQEEDLKKVMRRMQKEMDRGLRLETHEEASVKMLPTYVRSTPEGSEV
 GDFLSLDLGGTNFRVMLVKVGEEGQWSVKTKHQMYSIPEDAMTGAEMLFDYISECISDFLDKHQMKH
 KKLPLGFTFSFPVRHEDIDKGILLNWTGFKASGAEGNNVGLLRDAIKRRGDFEMDVVAMVNDTVATMI
 SCYYEDHQCEVGMIVGTGCNACYMEEMQNVELVEGDEGRMCVNTEWGAFGDSGELDEFLLEYDRLVDESS
 ANPGQQLYEKLIGGKYMGEVLRLVLLRLVDENLLFHGEASEQLRTRGAFETRFSQVESDTRKQIYNI
 LSTLGLRPSTTDCDIVRRACESVSTRAAHMCSAGLAGVINRMRESRSEDVMRITVGVDSVYKHLHPSFKE
 RFHASVRRLTPSCEITFIESEEGSGRGAALVSAVACKACMLGQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species:

Human

Serotype:

AAV-2

ACCN:

NM_033508

ORF Size:

1392 bp

Buffer:

PBS with 0.001% Pluronic F68

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033508.1</u>
RefSeq Size:	2566 bp
RefSeq ORF:	1395 bp
Locus ID:	2645
UniProt ID:	<u>P35557</u>
Cytogenetics:	7p13
MW:	51.9 kDa