

Product datasheet for **RC218941L3V**

TRPC4 (NM_016179) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRPC4 (NM_016179) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRPC4
Synonyms:	HTRP-4; HTRP4; TRP4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016179
ORF Size:	2931 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218941).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016179.1
RefSeq Size:	3440 bp
RefSeq ORF:	2934 bp
Locus ID:	7223
UniProt ID:	Q9UBN4
Cytogenetics:	13q13.3
Domains:	ANK, ion_trans
Protein Families:	Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane




[View online »](#)

MW: 111.9 kDa

Gene Summary: This gene encodes a member of the canonical subfamily of transient receptor potential cation channels. The encoded protein forms a non-selective calcium-permeable cation channel that is activated by Gq-coupled receptors and tyrosine kinases, and plays a role in multiple processes including endothelial permeability, vasodilation, neurotransmitter release and cell proliferation. Single nucleotide polymorphisms in this gene may be associated with generalized epilepsy with photosensitivity. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2011]

Product images:



[RC218941L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC218941L3V particle to overexpress human TRPC4-Myc-DDK fusion protein.