

Product datasheet for **RC218944L4V**

JHDM1D (KDM7A) (NM_030647) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	JHDM1D (KDM7A) (NM_030647) Human Tagged ORF Clone Lentiviral Particle
Symbol:	JHDM1D
Synonyms:	JHDM1D
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030647
ORF Size:	2823 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218944).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030647.1
RefSeq Size:	9178 bp
RefSeq ORF:	2826 bp
Locus ID:	80853
UniProt ID:	Q6ZMT4
Cytogenetics:	7q34
MW:	106.4 kDa



[View online »](#)

Gene Summary:

Histone demethylase required for brain development. Specifically demethylates dimethylated 'Lys-9' and 'Lys-27' (H3K9me2 and H3K27me2, respectively) of histone H3 and monomethylated histone H4 'Lys-20' residue (H4K20Me1), thereby playing a central role in histone code. Specifically binds trimethylated 'Lys-4' of histone H3 (H3K4me3), affecting histone demethylase specificity: in presence of H3K4me3, it has no demethylase activity toward H3K9me2, while it has high activity toward H3K27me2. Demethylates H3K9me2 in absence of H3K4me3. Has activity toward H4K20Me1 only when nucleosome is used as a substrate and when not histone octamer is used as substrate.[UniProtKB/Swiss-Prot Function]