

## Product datasheet for **RC218960L2V**

### TrkB (NTRK2) (NM\_001018065) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TrkB (NTRK2) (NM_001018065) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TrkB
Synonyms:	DEE58; EIEE58; GP145-TrkB; OBHD; trk-B; TRKB
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001018065
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218960).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001018065.1</a> , <a href="#">NP_001018075.1</a>
RefSeq Size:	8696 bp
RefSeq ORF:	1662 bp
Locus ID:	4915
UniProt ID:	<a href="#">Q16620</a>
Cytogenetics:	9q21.33
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	MAPK signaling pathway, Neurotrophin signaling pathway



[View online »](#)

**MW:** 60.99 kDa

**Gene Summary:** This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the MAPK pathway. Signalling through this kinase leads to cell differentiation. Mutations in this gene have been associated with obesity and mood disorders. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014]