

Product datasheet for **RC218963L3V**

C1S (NM_001734) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C1S (NM_001734) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1S
Synonyms:	EDSPD2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001734
ORF Size:	2064 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218963).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001734.2
RefSeq Size:	2714 bp
RefSeq ORF:	2067 bp
Locus ID:	716
UniProt ID:	P09871
Cytogenetics:	12p13.31
Domains:	CCP, CUB, Tryp_SPc, EGF_CA, EGF
Protein Families:	Druggable Genome, Protease


[View online »](#)

Protein Pathways: Complement and coagulation cascades, Systemic lupus erythematosus

MW: 76.68 kDa

Gene Summary: This gene encodes a serine protease, which is a major constituent of the human complement subcomponent C1. C1s associates with two other complement components C1r and C1q in order to yield the first component of the serum complement system. Defects in this gene are the cause of selective C1s deficiency. [provided by RefSeq, Mar 2009]