

## Product datasheet for **RC218966L2V**

### HIF-1 alpha (HIF1A) (NM\_181054) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HIF-1 alpha (HIF1A) (NM_181054) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HIF-1 alpha
Synonyms:	bHLHe78; HIF-1-alpha; HIF-1A; HIF-1alpha; HIF1; HIF1-ALPHA; MOP1; PASD8
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_181054
ORF Size:	2205 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218966).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_181054.2</a>
RefSeq Size:	3955 bp
RefSeq ORF:	2208 bp
Locus ID:	3091
UniProt ID:	<a href="#">Q16665</a>
Cytogenetics:	14q23.2
Protein Families:	Transcription Factors
Protein Pathways:	mTOR signaling pathway, Pathways in cancer, Renal cell carcinoma



[View online »](#)

MW: 83.2 kDa

**Gene Summary:** This gene encodes the alpha subunit of transcription factor hypoxia-inducible factor-1 (HIF-1), which is a heterodimer composed of an alpha and a beta subunit. HIF-1 functions as a master regulator of cellular and systemic homeostatic response to hypoxia by activating transcription of many genes, including those involved in energy metabolism, angiogenesis, apoptosis, and other genes whose protein products increase oxygen delivery or facilitate metabolic adaptation to hypoxia. HIF-1 thus plays an essential role in embryonic vascularization, tumor angiogenesis and pathophysiology of ischemic disease. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2011]